

For professionals

ReSound GN

# ReSound Vivia

Sound is more vivid than ever  
with Intelligence Augmented  
– a human approach to AI



GN Bringing People Closer

ReSound Vivia™ makes a world full of noise beautiful, by enhancing the natural hearing process and ability to focus on what matters. Let your clients embrace the vividness of life with the world's best solution for hearing in noise.<sup>1-4</sup>

# Intelligent solutions

## Intelligence Augmented

We believe that Artificial Intelligence can support the brain and its natural hearing abilities, but never replace it, as the brain is the most powerful sound processor. Our human approach to AI augments each person's unique strengths, so when extra help is needed, their brain can continue its natural sound processing.

## Intelligent Focus

Powered by our Deep Neural Network (DNN) chip\* – trained on real-life sounds and scenarios – this feature spotlights speech and contrasts unwanted noise for outstanding sound clarity and focus.

## Small, smaller, smallest



With our new advanced and efficient chip, you don't need to compromise on size when picking the best for your clients, such as the world's smallest AI hearing aid.<sup>4</sup>

## Built for wearing

ReSound Vivia gives a great wearing experience with all-weatherproof IP68-rated designs<sup>4</sup> and exchangeable mic filters for easy in-clinic servicing.

## All-day battery

Our microRIE style offers 30 hours of use on one charge – and 20 hours when streaming 50% of the time or using our new Intelligent Focus feature.<sup>4</sup>



## Future-proof connectivity

Offer the world's first Auracast Assistant for hearing aids<sup>4</sup> in the ReSound Smart 3D™ app to connect to Auracast™ broadcast audio, as well as our Auracast enabled TV-Streamer+ and Multi-Mic+.



1. Jespersen (2025)  
2. Cui & Groth (2024)  
3. Quilter (2021)  
4. GN Proprietary data on file

\* Available for Tech Level 9

# ReSound Vivia style options



Available in a beautiful selection of modern colours



## Key feature overview

For a detailed feature overview across styles, please refer to the data sheets

Features	VI9	VI7	VI5	VI4
<b>Sound Quality</b>				
WARP bands	★★★★ 17	★★★☆☆ 14	★★★☆☆ 12	★★★☆☆ 12
Environment recognition	★★★★ (Environmental Optimizer II)	★★★☆☆ (Environmental Optimizer)	★★★☆☆ (Environmental Optimizer)	★★★☆☆ (Environmental Optimizer)
Directionality	★★★★ (360 All-Around)	★★★☆☆ (All-Access)	★★★☆☆ (All-Access)	★★☆☆☆ (Binaural Directionality III)
Speech intelligibility in noise	★★★★ (Intelligent Focus / Clear Focus)	★★★☆☆ (Front Focus)	★★☆☆☆ (Ultra Focus)	★★☆☆☆ (AutoScope Adaptive Directionality)
Noise reduction	★★★★ (5 settings)	★★★☆☆ (4 settings)	★★★☆☆ (3 settings)	★★☆☆☆ (2 settings)
Impulse noise reduction	★★★★ (3 settings)	★★★★ (3 settings)	★★★☆☆ (2 settings)	★★☆☆☆ (2 settings)
Feedback management	★★★★ (DFS Ultra III)	★★★★ (DFS Ultra III)	★★★★ (DFS Ultra III)	★★★★ (DFS Ultra III)
Low frequency boost	★★★★ (3 settings)	★★★★ (3 settings)	★★★★ (3 settings)	★★★★ (3 settings)
<b>Connectivity</b>				
Hands-free calls	★★★★★ Hands-free calls for iPhone, iPad and compatible Android™ devices*			
Public connectivity	★★★★★ Auracast support for broadcast opportunities in private & public. Telecoil built into RIE 13 battery style, BTEs and in Multi-Mic+			
Accessories	★★★★★ TV-Streamer+, Multi-Mic+, Remote Control, ReSound Smart 3D app			
<b>Fitting Flexibility</b>				
Wireless CROS	★★★★★			
Gain handles	17	14	12	12
Programs	4			
Patient experience profile	Yes			
Flex Demo	Yes			
Bimodal Compatibility	ReSound Vivia is not bimodally compatible. Cochlear™ and GN continue to work together to bring state of the art bimodal solutions to Hearing Care Professionals and users under the Smart Hearing Alliance.			

\*Hands-free calls are compatible with iPhone 11 or later, iPad Pro 12.9-inch (5th generation), iPad Pro 11-inch (3rd generation), iPad Air (4th generation), iPad mini (6th generation), iPad (10th generation) or later, with software updates iOS 15.3 and iPadOS 15.3 or later. Bluetooth® Low Energy (LE) Audio works with ReSound Vivia and compatible mobile devices with Bluetooth 5.3 or later. See compatibility on [gnhearing.com/products/compatibility](https://gnhearing.com/products/compatibility)

**Worldwide Headquarters**  
GN Hearing A/S  
Lautrupbjerg 7  
DK-2750 Ballerup  
Denmark  
Tel.: +45 4575 1111  
[gnhearing.com](https://gnhearing.com)  
CVR no. 55082715

**United Kingdom**  
GN Hearing UK Ltd.  
Unit 13  
Talisman Business Centre  
Bicester  
Oxon OX26 6HR  
Tel.: +44 1869 352800  
[gnhearing.com](https://gnhearing.com)

**Australia**  
GN Hearing Australia Pty Ltd  
Suite 1001, Level 1  
1 Julius Avenue  
North Ryde, NSW 2113  
Australia  
Tel.: (free) 1800 658 955  
[gnhearing.com](https://gnhearing.com)

**New Zealand**  
GN Hearing New Zealand  
Limited  
Ground Floor, North Entrance  
4 Fred Thomas Drive  
Takapuna, Auckland, 0622  
Tel.: (free) 0800 900 126  
[gnhearing.com](https://gnhearing.com)